

# Dimension sheet torso

For orders please copy and send filled in form to:  
 E-Mail: info@albrechtgmbh.com  
 Fax: 08051-961 29 30

Order

Quotation

Date: \_\_\_\_\_

Medical supply store: \_\_\_\_\_

Customer ID: \_\_\_\_\_

Orthopaedic technician: \_\_\_\_\_

Patient ID: \_\_\_\_\_

## Product name

- MKS® Body belt
- MKS® Body belt plus
- MKS® Vario
- MKS® Lumbago
- MKS® Lumbago plus
- MKS® Lumbago Vario
- MKS® Pontsana®
- MKS® Pontsana® plus
- MKS® Pontsana® mobil
- MKS® Pontsana® stabil
- MKS® Osteo
- MKS® Osteo plus
- MKS® Osteo flex
- MKS® Osteo flex plus
- MKS® Lumbo
- MKS® Lumbo forte
- MKS® Lumbo plus
- MKS® Lumbo forte plus
- MKS® Lumbo mobil
- MKS® Lumbo fenestra
- MKS® Extensa
- MKS® Extensa plus
- MKS® Thorecta
- MKS® Thorecta plus
- MKS® Core Compression Bandage
- other: \_\_\_\_\_

**Abdominal height**  
 Maximum front height of the textile body belt

**Front height**  
 Measured from the upper edge of the pubic bone to the sternum (3-4 below the collarbone)

**Body circumference**  
 Measured at the middle of the torso (just above the navel)

**Back height reclinatio**  
 Measured from the sacrum (beginning of gluteal fold) to circa 2-3 cm below cervical vertebra C7

**Back height**  
 Measured from the sacrum to circa 2-3 cm below the tips of the shoulder blades



### Technical implementation (to be filled in by albrecht)

- Size modification
- Modular adaptation

Art.-No.: \_\_\_\_\_

\_\_\_\_\_

Approval: \_\_\_\_\_